





LINA is a compact relaxation wingback lounge chair, perfect for long-term sitting. In addition to the fixed lounge chair version, it is also available in manual and electric reclining options. The transfer version is also available.

Use

Suitable for relaxation areas.

Seat and backrest

Internal structure made of solid wood, with a seat and backrest suspension system using highly durable elastic straps. Covered with flexible, high-density, and HR polyurethane foam. Continuous upholstery on the seat and footrest to prevent gaps that could cause injuries.

Armrests

Internal structure made of solid wood covered with flexible, high-density, and HR polyurethane foam. Optionally, the upholstered arms can be raised to facilitate lateral transfer of the user.

Legs

Made of solid wood, varnished in any of the colours from our sample collection or in any colour specified by the client.

Mechanisms

Manual recliner chair: the user can manually adjust the coordinated inclination of the backrest and footrest. Optionally, it can include a lever to facilitate the opening of the armchair.

Manual recliner chair + ZERO GRAVITY position: the user can manually adjust the coordinated inclination of the backrest and footrest.

It includes the possibility of adopting the GRAVITY ZERO position. In the GRAVITY ZERO position, the legs are at the same level as the heart to promote blood circulation. It is a relaxing position where the lumbar and cervical muscles relax, as less pressure is exerted on them.

Electric recliner chair: the user can adjust the desired inclination thanks to a remote control. Electrically powered.

During the opening, by keeping the button pressed, you can complete the movement to the final relaxation position or interrupt it at any intermediate points. Once closed, it returns to the sitting position.

Electric recliner chair + lift – ZERO GRAVITY position: the user can adjust the desired inclination using a remote control. Electrically powered.

During the opening, by keeping the button pressed, you can complete the movement to the final relaxation position or interrupt it at any intermediate points. Once closed, using the same button, the user can continue the movement until the standing position.

It includes the possibility of adopting the GRAVITY ZERO position. In the GRAVITY ZERO position, the legs are at the same level as the heart to promote blood circulation. It is a relaxing position where the lumbar and cervical muscles relax, as less pressure is exerted on them.

Transfer

The TRANSFER versions incorporate a transfer handle and high-quality 10mm diameter wheels, two of them with brakes, allowing the chair to be locked when not in use.

Upholstery

Choose from our wide range of waterproof vinyl fabrics with antimicrobial, antibacterial, and antifungal protection.

Certifications

Seniorcare Interiors holds the ISO 9001 and ISO 14001 certifications, attesting that its work processes maintain a high standard of quality, covering the entire product lifecycle from the design phase to production and subsequent distribution of the product.

Please consult the SENIORCARE finishing chart for information about the certifications of the fabrics.



Seniorcare Interiors is certified to ISO 9001 and ISO 14001 standards.

LINA manual recliner armchair

CODE

height	width	depth	depth abierto	seat height	arm. height
112	70	82	170	46	62

LIN0011TA



* Fabric yardage: 4,7 m

LINA manual recliner ZERO GRAVITY

CODE

height	width	depth	depth abierto	seat height	arm. height
112	70	82	170	46	62

LIN0020TA



* Fabric yardage: 4,7 m

LINA fixed armchair

CODE

height	width	depth	depth abierto	seat height	arm. height
112	70	82	170	46	62

LIN0002TA



* Fabric yardage: 4,7 m

LINA electric recliner

CODE

height	width	depth	depth abierto	seat height	arm. height
112	70	82	170	46	62

LIN0013TA



* Fabric yardage: 4,7 m

LINA electric lift recliner ZERO GRAVITY

CODE

height	width	depth	depth abierto	seat height	arm. height
112	70	82	170	46	62

LIN0030TA



* Fabric yardage: 4,7 m

Increase for transfer option

CODE

LIN0110TA

